

Архангельск (8182)63-90-72
Астана (7172)727-132
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13

Сургут (3462)77-98-35
Тверь (4822)63-31-35
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Ярославль (4852)69-52-93

www.ammed.nt-rt.ru || abm@nt-rt.ru


МАНЕКЕНЫ И НАБОРЫ ЛИЦЕВЫХ ПЛАСТИН

Техническое описание

® CPR Pal

SPARE PARTS

	CAT. NO.
1 ACD CPR base board	259 000 701
2 Set of 100 head bags	174 000 701
3 Set of 25 mouth/nose pieces	071 000 586
4 Carrying bag for CPR Pal	259 000 703
5 Locking cotters of 4 pcs. and 2 screws	268 000 513
6 Spring	268 000 520



Key Features

- Wireless connections
- Ambu's unique hygienic system
- The rigidity of the chest is adjustable
- Airways open only when head is correctly extended
- Jaw thrust possible
- PC compatible
- Training mat integrated in the carrying bag

® Man W + I

The Man is a highly developed instruction- and training manikin for the simulation of realistic and correct anatomical and physiological conditions during cardiopulmonary resuscitation by means of rescue breathing and chest compression.

Realistic scenarios

The Man W (Wireless) does not need cables and wires. Which means it can be placed in realistic surroundings away from computers and power supplies. This gives a great amount of freedom during training and because there are no wires or cables lying around it gives an even more realistic scenario.

Low risk of cross-contamination

The Man incorporates an unique, patented hygienic system under which all trainees get their own face piece and head bag, thereby minimizing any risk of cross-infection between the trainees. The hygienic system also features realistic expiration from the Man through the mouth and nose.

Mechanical Instrument

Version W:

The Man W models are equipped with a mechanical monitoring instrument as well as digital sensors showing effectiveness of the resuscitation. It also shows results achieved for insufflations volume and depth of external chest compression, as well as it indicates any stomach inflation and incorrect hand positioning.

Man W models can be connected to a computer through USB or WLAN using the CPR-Software for a more detailed analysis of the performed resuscitation.

Version I:

The Man I (Instrument) models are equipped with a mechanical monitoring instrument showing effectiveness of the resuscitation, including results achieved for ventilation volume and depth of external chest compression. Moreover it indicates any stomach inflation and incorrect hand positioning.

® SAM

Description

The SAM is an inexpensive training manikin with all the well-known and sought-after look-listen-feel features. The SAM incorporates the unique and patented hygienic system that enables training with a face piece and head bag. Especially the closed shape of the torso with its anatomical correct landmarks enables students to learn the complete algorithm of Basic Life support including the use of Automatic External Defibrillators (AED) and the placement of the AED Pads.

Open Airway

The airway is opened only when the head is tilted - there are realistic airway resistance and chest rise.

Hygienic System

SAM has got the same unique hygienic solution as the Man with disposable head bags and washable face pieces.

Compact Bag

The head of the Sam can be placed inside the manikins so everything is packed in only one small bag.



ACCESSORIES

Head bags (sets of 100)	234 000 702
Face piece (sets of 5)	234 000 703

SPECIFICATIONS

Weight:	2 kg
Length:	70 cm
Width	33 cm
Height	20 cm
Max. insuf. volume	1.8 liters
Max compression depth	7,5 cm
Storage temperature	-20° to +70° C




Compact storage

ORDERING INFORMATION

SAM with 25 head bags and 1 face piece, supplied in a soft carrying bag	A268 003 000
Head bag (supplied in sets of 100)	234 000 702
Face piece (supplied in sets of 5)	234 000 703



Compression depth indicator



Placemend of AED pads

ACCESSORIES

5 Face pieces	234 000 703
100 Head bags	234 000 702
Arms complete	234 000 704
Arms, legs, jogging suit trousers and bag	234 000 706
Jogging suit jacket	234 000 790
Jogging suit trousers	234 000 792
Legs with trousers	234 000 786
Carrying bag	234 000 794

SPECIFICATIONS

Weight: Torso with carrying case	12 kg
Full body with carrying case	17 kg
Total length: Torso	80 cm
Full body	170 cm

ORDERING INFORMATION Man Model I,

Torso	A234 204 000
Man Model W, Torso with CPR-Software	A234 307 000
Upgrade Kit - Full Body	234 000 706
Upgrade Kit - Arms	234 000 704

MATERIALS

Basic unit	Polyethylene
Instrument part	ABS plastic
Skull	PVC, hard
Chest, face	PVC, soft



Man computer interface



Arms and legs for Man



Bag for Man

Архангельск (8182)63-90-72
Астана (7172)727-132
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Волгоград (844)278-03-48
Волгодга (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13

Сургут (3462)77-98-35
Тверь (4822)63-31-35
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Ярославль (4852)69-52-93